


CODE	OUTSIDE	INSIDE	WEIGHT
11-700	595 X 365 X 125 mm	530 X 350 X 100 mm	2.620 kg
11-701	595 X 365 X 187 mm	530 X 350 X 162 mm	3.440 kg



CODE	OUTSIDE	INSIDE	WEIGHT
11-702	595 X 365 X 105 mm	530 X 350 X 100 mm	2.20 kg
11-703	595 X 365 X 167 mm	530 X 350 X 167 mm	3.00 kg




Lid Barrier system




11-771	
11-772	
11-773	
11-774	
11-775	
11-776	




GREY
YELLOW
GREEN
BLUE
RED
BLACK




1. Before sterilization the valves are in relax position.




2. During the sterilization by effect of vacuum the valves are opened completely and the dry air goes out of the container.



3. At the end the sterilization cycle the valves are opened in the reverse way by effect of vacuum and the sterilization cycle is completed.




4. At the end of the sterilization cycle the sterilizer equates the pressure level and the valves of the Barrier Model Container are closed.



5. After sterilization the valves are in relax position.



Sterile Container System



Dental Mini & Special Model

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)
320 X 190 X 50 mm 	10-1800  10-1801  10-1802  10-1803  10-1804  10-1805 	10-1820  10-1821  10-1822  10-1823  10-1824  10-1825 	10-1840  10-1841  10-1842  10-1843  10-1844  10-1845 	10-1860 	10-1865 
320 X 190 X 65 mm 	10-1806  10-1807  10-1808  10-1809  10-1810  10-1811 	10-1826  10-1827  10-1828  10-1829  10-1830  10-1831 	10-1846  10-1847  10-1848  10-1849  10-1850  10-1851 	10-1861 	10-1866 
320 X 190 X 130 mm 	10-1812  10-1813  10-1814  10-1815  10-1816  10-1817 	10-1832  10-1833  10-1834  10-1835  10-1836  10-1837 	10-1852  10-1853  10-1854  10-1855  10-1856  10-1857 	10-1862 	10-1867 

Ps: All container measurements cover the outer size of the lids.



13-1070 274 x 172 x 43 mm
WIRE BASKET



13-611 274 x 172 x 43 mm
WIRE BASKET WITH LID.



13-609 273 x 180 x 13 mm
STORAGE FRAME
for dental model




12-529 260 x 180 mm
SILICONE MAT



13-610 277 x 175 mm
DRAPE RETAINER

Perforated Bottoms




10-1860

Solid Bottoms




10-1865



10-1861



10-1866



10-1862



10-1867

Perforated Lids



13-615


Single use paper filter 100 pcs.
With indicator No marking.




12-743

Teflon filter
Number of sterilization cycle > 1200

Bio-Barrier Lids



Removal Lids



Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)
315 X 135 X 40 mm 	10-1900  10-1901  10-1902  10-1903  10-1904  10-1905 	10-1920  10-1921  10-1922  10-1923  10-1924  10-1925 	10-1940  10-1941  10-1942  10-1943  10-1944  10-1945 		
315 X 135 X 70 mm 	10-1906  10-1907  10-1908  10-1909  10-1910  10-1911 	10-1926  10-1927  10-1928  10-1929  10-1930  10-1931 	10-1946  10-1947  10-1948  10-1949  10-1950  10-1951 	10-1960 	10-1965 
315 X 135 X 100 mm 	10-1912  10-1913  10-1914  10-1915  10-1916  10-1917 	10-1932  10-1933  10-1934  10-1935  10-1936  10-1937 	10-1952  10-1953  10-1954  10-1955  10-1956  10-1957 	10-1961 	10-1966 
				10-1962 	10-1967 

Ps: All container measurements cover
the outer size of the lids.



13-1073 267 x 125 x 50 mm
WIRE BASKET



12-993 267 x 125 x 50 mm
WIRE BASKET WITH LID




12-990 270 x 125 mm
DRAPE RETAINER



12-991
Portable carrier for small containers.
Wire suspenders adjustable as required.
Adjusting key of suspender is included.
Width: 325 mm, Height: 328 mm, Depth: 145 mm



12-992 270 x 125 mm
WIRE BASE



12-520 275 x 125 mm
SILICONE MAT

Perforated Bottoms



10-1960



10-1961



10-1962

Solid Bottoms



10-1965



10-1966



10-1967



12-761

Single use paper filter 100 pcs.
With indicator. No marking.



12-744

Teflon filter
Number of sterilization scycle > 1200

Perforated Lid



Removal Lids



IMPLANT Model Containers

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)	NON FILTER
	 Box (Non-Perforated)	 Box (Perforated)	 Box (Non-Perforated)	
300 X 149 X 65 mm 	09-010  09-011  09-012  09-013  09-014  09-015 	09-022  09-023  09-024  09-025  09-026  09-027 	09-034  09-035  09-036  09-037  09-038  09-039 	
500 X 169 X 75 mm 	09-016  09-017  09-018  09-019  09-020  09-021 	09-028  09-029  09-030  09-031  09-032  09-033 	09-040  09-041  09-042  09-043  09-044  09-045 	

Ps: All container measurements cover the outer size of the lids.

MINI IMPLANT Model Containers

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)	NON FILTER
	 Box (Non-Perforated)	 Box (Perforated)	 Box (Non-Perforated)	
160 X 70 X 60 mm 	09-200  09-201  09-202  09-203  09-204  09-205 	09-212  09-213  09-214  09-215  09-216  09-217 	09-224  09-225  09-226  09-227  09-228  09-229 	
200 X 145 X 60 mm 	09-206  09-207  09-208  09-209  09-210  09-211 	09-218  09-219  09-220  09-221  09-222  09-223 	09-230  09-231  09-232  09-232  09-234  09-235 	

Ps: All container measurements cover the outer size of the lids.

**BRUSH Model
Containers**


-  **12-044** Black
-  **12-045** Red
-  **12-046** Blue
-  **12-047** Green
-  **12-048** Yellow
-  **12-049** Grey

Container for brush 280 x 285 x 85 mm with wall mounting fixtures, Maximum 20 each brushes can be stored container to take the brush, please move the handle to right or left direction under the container.

non filter

Specifications / Sizes	Lid (Perforated)	Lid (Perforated)
	Box (Non-Perforated)	Box (Perforated)
250 X 60 X 50 mm 	09-060 09-061 09-062 09-063 09-064 09-065	09-072 09-073 09-074 09-075 09-076 09-077
450 X 70 X 70 mm 	09-066 09-067 09-068 09-069 09-070 09-071	09-078 09-079 09-080 09-081 09-082 09-083

Ps: All container measurements cover the outer size of the lids.



SCOPY Model Containers

Specifications / Sizes	Lid (Perforated)	Lid (Perforated)
	Box (Non-Perforated)	Box (Perforated)
250 X 60 X 50 mm 	09-100 09-101 09-102 09-103 09-104 09-105	09-115 09-116 09-117 09-118 09-119 09-120
450 X 70 X 70 mm 	09-106 09-107 09-108 09-109 09-110 09-111	09-121 09-122 09-123 09-124 09-125 09-126

Ps: All container measurements cover the outer size of the lids.



SCOPY Model Containers
(with silicon stand, lock, label)

Specifications Sizes	Lid (Perforated)	Lid (Perforated)	Lid (Non Perforated)		
	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)	Box (Perforated)	Box (Non-Perforated)
285 X 280 X 55 mm 	10-1710  10-1711  10-1712  10-1713  10-1714  10-1715 	10-1730  10-1731  10-1732  10-1733  10-1734  10-1735 	10-1750  10-1751  10-1752  10-1753  10-1754  10-1755 		
285 X 280 X 85 mm 	10-1716  10-1717  10-1718  10-1719  10-1720  10-1721 	10-1736  10-1737  10-1738  10-1739  10-1740  10-1741 	10-1756  10-1757  10-1758  10-1759  10-1760  10-1761 	10-1770 	10-1775 
				10-1771 	10-1776 

Ps: All container measurements cover
the outer size of the lids.




12-490 255 x 250 x 30 mm
12-491 255 x 250 x 50 mm
 WIRE BASKET




12-482 255 x 255 x 30 mm
 WIRE BASE



12-480 255 x 255 mm
 DRAPE RETAINER



12-481 255 x 255 x 30 mm
 DRAPE RETAINER



12-521 255 x 255 mm
 SILICON MAT




12-740
 Single use paper filter 100 pcs.



12-742
 Teflon Filter
 Number of sterilization cycle > 1200

Perforated Bottoms



10-1770

Solid Bottoms



10-1775



10-1771



10-1776


Perforated Lids



Bio-Barrier Lids



Safety Lids



Bio-Barrier Safety Lids



Removal Lids



Container System & Baskets

Classic & Bio-Barrier Model



 **Thempson**
Surgical Instruments

Surgical innovation Sytems
From Access to Repair

Address:

Veko Giz Plaza, Meydan Sokak No.3/43, Kat.13 Maslak, Istanbul,Turkey

Phone: +90 212 705 36 16, Fax: +90 212 705 36 36

info@thempson.com

www.thempson.com

Thempson[®] is registered trademark of Sueda Danışmanlık ve Dış Tic. Ltd. Şti. All Rights are Reserved.